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

Parameters command (Setup menu)

This command opens a dialog box to set experimental parameters for the technique currently selected in the [Electrochemical Techniques dialog box](#).

This command has a toolbar button:

The following parameter dialog boxes available for the CHI600E series are listed below. For details about each technique, please refer to the description of each dialog box.


- Parameters for Cyclic Voltammetry
- Parameters for Linear Sweep Voltammetry
- Parameters for Staircase Voltammetry
- Parameters for Tafel Plot
- Parameters for Chronoamperometry
- Parameters for Chronocoulometry
- Parameters for Differential Pulse Voltammetry
- Parameters for Normal Pulse Voltammetry
- Parameters for Differential Normal Pulse Voltammetry
- Parameters for Square Wave Voltammetry
- Parameters for A.C. Voltammetry
- Parameters for 2nd Harmonic A.C. Voltammetry
- Parameters for FT A.C. Voltammetry
- Parameters for Amperometric i-t Curve
- Parameters for Differential Pulse Amperometry
- Parameters for Double Differential Pulse Amperometry
- Parameters for Triple Pulse Amperometry
- Parameters for Integrated Pulse Amperometric Detection
- Parameters for Bulk Electrolysis with Coulometry
- Parameters for Hydrodynamic Modulation Voltammetry
- Parameters for Sweep-Step Functions
- Parameters for Multi-Potential Steps
- Parameters for A.C. Impedance
- Parameters for Impedance - Time
- Parameters for Impedance - Potential
- Parameters for Chronopotentiometry
- Parameters for Chronopotentiometry with Current Ramp
- Parameters for Multi-Current Steps
- Parameters for Potentiometric Stripping Analysis
- Parameters for Electrochemical Noise Measurement
- Parameters for Open Circuit Potential - Time

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The screenshot shows the 'System Setup' dialog box with the following settings and annotations:

- Communication Port:** A list of ports from Com 1 to Com 12. **Annotation:** A red bracket on the left indicates '通讯端口应与电脑设备管理器端口显示号码一致' (Communication port should be consistent with the port number displayed in the computer's device manager).
- Line Frequency:** Radio buttons for 60 Hz and 50 Hz. **Annotation:** '电源频率 50Hz' (Power frequency 50Hz) with an arrow pointing to the 50 Hz option.
- Potential Axis:** Radio buttons for Positive Left and Positive Right. **Annotation:** A red bracket on the right indicates '坐标轴方向: 正值向右, 正值向上' (Coordinate axis direction: positive value to the right, positive value up).
- Current Axis:** Radio buttons for Positive Up and Positive Down.
- Com Port Speed:** Radio buttons for Standard and Fast.
- Current Polarity:** Radio buttons for Cathodic Positive and Anodic Positive. **Annotation:** A blue arrow points to the Anodic Positive option with the text '电流极性, 一般为阳极电流为正' (Current polarity, generally anodic current is positive).
- Windows:** Radio buttons for English and Oriental.
- Data Length:** A dropdown menu set to 8192K.
- Checkboxes:** 'Save retrieve data during run', 'Present Data Override Warning', 'Save Text File As Well', and 'Check COM port while idling'. **Annotation:** A red bracket on the right groups these checkboxes with the text '运行中保存数据(断电备份) 数据覆盖提示 同时保存txt数据' (Save data during run (power off backup) Data override warning Also save txt data).
- Buttons:** OK, Cancel, Help, Set Com Port, and Erase ADC Calibration Coefficients.

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Quiet Time (sec)

Sensitivity (A/V)

☐ Open Circuit Potential as the Center Potential

☐ Auto Sens if Scan Rate ≤ 0.01 V/s

☐ Auxiliary Signal Recording if Scan Rate ≤ 0.25 V/s

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


Parameter	Range	Description
Init E (V)	-10 - +10	Initial potential
Final E (V)	-10 - +10	Final potential
Incr E (V)	1e-3 - 0.05	Increment potential of each step
Sweep Segments	1 - 2	Sweep segments; each segment is half cycle
Scan Rate (V/s)	1e-6 - 0.01	Potential scan rate
Quiet Time (sec)	0 - 100000	Quiescent time before potential scan
Sensitivity (A/V)	1e-12 - 0.1	Sensitivity scale
Auto Sens	Check or Uncheck	Automatic sensitivity switching during run

Chronoamperometry (计时电流法)

Chronoamperometry Parameters

Init E (V)

High E (V)

Low E (V)

Initial Step Polarity.....

Number of Steps

Pulse Width (sec)

Sample Interval (sec)

Quiet Time (sec)

Sensitivity (A/V)

☐ Auxiliary Signal Recording When Sample Interval >= 0.005 s

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



Differential Normal Pulse Voltammetry

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Differential Normal Pulse Voltammetry

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



Parameter	Range	Description
Init E (V)	-10 - +10	Initial potential
Sample Interval (s)	1e-6 - 50	Data sampling interval
Run Time (s)	0.001 - 500000	Total measurement time
Quiet Time (s)	0 - 100000	Quiescent time before potential step
Scales During Run	1, 2, 3	Number of current display scales
Sensitivity (A/V)	1e-12 - 0.1	Sensitivity scale
Aux. Signal Rec.	Check or Uncheck	Record external signal when sample

Differential Pulse Amperometry (差分脉冲电流法)

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
Parameter	Range	Description
Init E (V)	-10 - +10	Initial potential
Cleaning E (V)	-10 - +10	Electrode cleaning potential
Cleaning Time (s)	0 - 32	Electrode cleaning time
Pulse E1 (V)	-10 - +10	First pulse potential
Pulse T1 (s)	0.01 - 32	First pulse time
Pulse E2 (V)	-10 - +10	Second pulse potential
Pulse T2 (s)	0.01 - 32	Second pulse time
Number of Cycles	10 - 100000	Number of Repetitive Cycles
Quiet Time (s)	0 - 100000	Quiescent time before potential pulses begin
Scales During Run	1, 2, 3	Number of current display scales
Sensitivity (A/V)	1e-12 - 0.1	Sensitivity scale
Open Circ. Clean.	Check or Uncheck	Cleaning step held at constant or open circuit potential

Double Differential Pulse Amperometry (双差分脉冲电流法)

Double Differential Pulse Amperometry

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
Triple Pulse Amperometry (三脉冲电流法)

Triple Pulse Amperometry Parameters

Potential 1
 E (V) 0
 Duration (s) 0.1
☐ Open Circuit

Potential 2
 E (V) 0
 Duration (s) 0.1


Potential 3
 E (V) 0
 Duration (s) 0.1
 Incr E (V) 0

Init E (V) 0
 Final E (V) 0
 Number of Cycles 400
 Quiet Time (sec) 0
 Scales during Run 1
 Sensitivity (A/V) 1.e-006

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Parameter	Range	Description
E1 (V)	-10 - +10	First pulse potential
Duration 1 (s)	0 - 32	First pulse duration
Open Circuit	Check or Uncheck	Hold step 1 at open circuit or constant potential
E2 (V)	-10 - +10	Second pulse potential
Duration 2 (s)	0 - 32	Second pulse duration
E3 (V)	-10 - +10	Third pulse potential
Duration 3 (s)	0.01 - 32	Third pulse duration
Incr E (V)	0 - 0.02	Increment potential for third step
Init E (V)	-10 - +10	Initial potential during quiescent time
Final E (V)	-10 - +10	Final potential for scan
Number of Cycles	10 - 100000	Number of Repetitive Cycles
Quiet Time (sec)	0 - 100000	Quiescent time before taking data
Scales During Run	1, 2, 3	Number of current display scales
Sensitivity (A/V)	1e-12 - 0.1	Sensitivity scale


Integrated Pulsed Amperometric Detection

Integrated Pulsed Amperometric Detection Parameters

Step 1: Start Start E (V) <input type="text" value="0"/> Hold Time (s) <input type="text" value="0.2"/>	Step 4: Hold Hold E (V) <input type="text" value="0"/> Hold Time (s) <input type="text" value="0.2"/>	OK Cancel Help
Step 2: Forward Scan Peak E (V) <input type="text" value="0"/> Scan Time (s) <input type="text" value="0.2"/>	Step 5: Oxidation Oxd E (V) <input type="text" value="0"/> Oxd Time (s) <input type="text" value="0.2"/>	
Step 3: Reverse Scan Return E (V) <input type="text" value="0"/> Scan Time (s) <input type="text" value="0.2"/>	Step 6: Reduction Red E (V) <input type="text" value="0"/> Red Time (s) <input type="text" value="0.2"/>	

No. of Cycles
 Quiet Time (sec)
 Scales during Run ...
 Sensitivity (A/V)

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Bulk Electrolysis with Coulometry (控制电位电解库仑法)

Bulk Electrolysis with Coulometry Parame... [X]

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
Measurement Mode above 100 Hz

☐ FT ☒ Single Freq. ☐ Galvanostatic


Measuring Time or Cycles and Points

Freq Range	Avg/Cycles	Points / Decade Freq
100K - 1M Hz :	<input type="text" value="1"/>	<input type="text" value="12"/>
10K - 100K Hz :	<input type="text" value="1"/>	<input type="text" value="12"/>
1K - 10K Hz :	<input type="text" value="1"/>	<input type="text" value="12"/>
100 - 1K Hz :	<input type="text" value="1"/>	<input type="text" value="12"/>
10 - 100 Hz :	<input type="text" value="1"/>	<input type="text" value="12"/>
1 - 10 Hz :	<input type="text" value="1"/>	<input type="text" value="12"/>
0.1 - 1 Hz (cycles)	<input type="text" value="1"/>	<input type="text" value="12"/>
0.01 - 0.1 Hz (cycles)	<input type="text" value="1"/>	<input type="text" value="12"/>
0.001 - 0.01 Hz (cycles)	<input type="text" value="1"/>	<input type="text" value="12"/>
.0001 - .001 Hz (cycles)	<input type="text" value="1"/>	<input type="text" value="12"/>
.00001 - .0001 Hz (cyc)	<input type="text" value="1"/>	<input type="text" value="12"/>

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
Parameter	Range	Description
Init E (V)	-10 - +10	Initial potential
Amplitude (V)	0.0001 - 1.5	A.C. amplitude
Frequency (Hz)	0.0001 - 1e6	A.C. frequency
Sample Interval (s)	1 - 20000	Data sampling interval
Run Time (s)	1 - 500000	Total experiment time
Cycles (< 10 Hz)	1 - 100	Number of repetitive cycles at each frequency
Quiet Time (s)	0 - 100000	Quiescent time before sampling data
Bias DC Current	off - range - on	Enable DC current bias during run
Sensitivity Scale	1e-12 - 0.1	Sensitivity scale (automatic or manual)

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
Data Storage Intvl (sec)

Number of Segments


Current Switching Priority
☒ Potential ☐ Time

Auxiliary Signal Recording When Sample Interval ≥ 0.0005 s
☒ No ☐ Ext Voltage ☐ Current Control

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Parameter	Range	Description
Cathodic Curr. (A)	-2 - +2	Controlled cathodic current
Anodic Current (A)	-2 - +2	Controlled anodic current
High E Limit (V)	-10 - +10	High potential limit
High E Hold T. (s)	0 - 100000	High potential hold time
Low E Limit (V)	-10 - +10	Low potential limit
Low E Hold T. (s)	0 - 100000	Low potential hold time
Cathodic Time (s)	0.005 - 100000	Cathodic run time
Anodic Time (s)	0.005 - 100000	Anodic run time
Initial Polarity	Cathodic/Anodic	Polarity for the first segment
Data Stor. Intvl (s)	0.0001 - 500	Data storage interval
Num. of Segments	1 - 1000000	Number of half cycles
Curr. Switch. Prior.	Potential or Time	Current polarity switching control
Aux. Signal Rec.	Check or Uncheck	Simult. record external signal when samp.

Chronopotentiometry with Current Ramp (电流扫描计时电位法)

Chronopotentiometry with Current Ramp... ✕

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


Potentiometric Stripping Analysis (电位溶出分析)

Potentiometric Stripping Analysis Para... ✕

Deposition E (V)	<input type="text" value="0"/>	<input type="button" value="OK"/> <input type="button" value="Cancel"/> <input type="button" value="Help"/>
Deposition Time (sec)	<input type="text" value="15"/>	
Final E (V)	<input type="text" value="0"/>	
Stripping Current (A)	<input type="text" value="0"/>	
Sample Interval (sec)	<input type="text" value="0.01"/>	
Quiet Time (sec)	<input type="text" value="0"/>	

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Parameter	Range	Description
Deposition E (V)	-10 - +10	Deposition potential
Deposit. Time (s)	0 - 1000000	Deposition time
Final E (V)	-10 - +10	Final potential (see Notes)
Stripping Curr. (A)	0 - 2	Controlled stripping current
Sample Interval (s)	0.0001 - 50	Sampling interval
Quiet Time (s)	0 - 100000	Quiescent time before taking data

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


Control Level of Stir Line
☒ Active High ☐ Active Low

Immediate Stir
Stir Time (sec) 10
Stir

Immediate Purge
Purge Time (sec) 240
Purge

Immediate Cell On
Potential (V) 0
Time (sec) 10
Cell On

SMDE Drops Before Run 1

Stabilizing Capacitor
☒ Automatic Setting
☐ On ☒ Off

Mercury Drop Collection
Number of Drops 50
Time Interval (sec) 1
Collect

OK
Cancel
Help

☐ 4 Electrodes
☒ Stir Between Run
☒ Purge Between Run
☒ Stir and Purge Simultaneously
☐ Test with Internal Dummy Cell
☐ Cell On Between Runs
☒ Return to Init E after Run

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


The 'Step Functions' dialog box contains the following fields and buttons:

- Start E (V) 0
- Duration (s) 0
- Step E1 (V) / i1 (A) ... 0
- Step Time 1 (s) 0.1
- Step E2 (V) / i2 (A) 0
- Step Time 2 (s) 0.1
- Step Segment 2
- ☐ Galvanostatic Mode
- Buttons: OK, Cancel, Help, Run

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

Parameter	Range	Description
Start E (V)	-10 - +10	Starting potential
Duration (s)	0 - 100000	Duration at Start E
Step E1 (V) / i1 (A)	-10 - +10	First step potential
Step Time 1 (s)	.0001 - 100000	Duration of each step
Step E2 (V) / i2 (A)	-10 - +10	Second step potential
Step Time 2 (s)	.0001 - 100000	Duration of each step
Step Segments	1 - 2000000000	Step segments; each segment is a half cycle
Galvanostatic Mode	Check or Uncheck	Galvanostatic or potentiostatic (unchecked) mode

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



Peak Definition


Peak Shape	Report Option
<input type="radio"/> Gaussian	<input checked="" type="checkbox"/> Peak or wave potential
<input checked="" type="radio"/> Diffusive	<input type="checkbox"/> Half peak potential
<input type="radio"/> Sigmoidal	<input checked="" type="checkbox"/> Peak or wave current
	<input checked="" type="checkbox"/> Peak area

Peak or wave search potential range (V) 0.5

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Calibration Curve dialog box. It contains a table for inputting standard concentrations and peak heights, fields for slope, intercept, and correlation, and axes title/unit fields.

	Concentration	Peak Height
Standard 1:	0	0
Standard 2:	0	0
Standard 3:	0	0
Standard 4:	0	0
Standard 5:	0	0
Standard 6:	0	0
Unknown:	0	0

Slope: 0
Intercept: 0
Correlation: 0

X Axis Title:
X Axis Unit:
Y Axis Title:
Y Axis Unit:
Header:
Note:

Buttons: OK, Cancel, Help, Read, Save, Calculate, Plot


Standard Addition dialog box. It contains a table for inputting unknown peak height and standard additions, fields for slope and correlation, and axes title/unit fields.

	Concentration	Peak Height
Unknown:	0	0
Addition 1:	0	0
Addition 2:	0	0
Addition 3:	0	0
Addition 4:	0	0

Slope: 0
Correlation: 0

X Axis Title:
X Axis Unit:
Y Axis Title:
Y Axis Unit:
Header:
Note:

Buttons: OK, Cancel, Help, Read, Save, Calculate, Plot

Data File Report [X]

	Species 1	Species 2	Species 3	Species 4
Species:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Ep From (V):	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>
Ep To (V):	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>
Slope:	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>
Intercept:	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>
Peak Shape:	<input type="text" value="Gaussian"/>			
Data Type:	<input type="text" value="Original"/>			
Number of Species:	<input type="text" value="1"/>			
Concentration Unit:	<input type="text"/>			
Report Filename:	<input type="text"/>			
Header:	<input type="text"/>			
Note:	<input type="text"/>			

Report Data Type
☒ Conc. ☐ Current

Buttons: OK, Cancel, Help, Read, Save, Report

Time Dependence [X]

Ep From (V):
Ep To (V):
Slope:
Intercept:
Peak Shape:
Data Type:
Report Name:
Y Axis Title: Y Axis Unit:
Header:
Note:

Report Data Type
☒ Conc. ☐ Current


Buttons: OK, Cancel, Help, Read, Save, Report, Plot

Combine Multiple Time based Data Files into One Time Based File

Folder and Base Filename:
Number of Files: Folder Browse
☒ Combine to One Point per File
Combine

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Cancel

Help

Cathodic Tafel Slope Potential Range

-0.389 to -0.329

Anodic Tafel Slope Potential Range

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


Data Proc Performed:

- Smoothing
- 1st Derivative
- 2nd Derivative
- 3rd Derivative
- 4th Derivative
- 5th Derivative
- Integration
- Semi-Derivative
- Semi-Integral
- Interpolation
- Baseline Correction
- Data Point Removing
- Data Point Modifying
- Bkgnd Subtraction
- Signal Averaging
- X Math Operation
- Y Math Operation

Header:

Note:

OK
Help

Data List


Jan. 25, 2016 22:12:10
Cyclic Voltammetry
File: baseline corr unsaved
Data Source: Experiment
Instrument Model: CHI760E
Header:
Note:

Init E (V) = -0.7
High E (V) = 0.3
Low E (V) = -0.7
Init P/N = P
Scan Rate (V/s) = 0.1
Segment = 474
Sample Interval (V) = 0.01
Quiet Time (sec) = 2
Auto Sens

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

- AC Impedance (IMP)
- Impedance - Time (IMPT)
- Impedance - Potential (IMPE)
- AC (including phase-selective) Voltammetry (ACV)
- Second Harmonic AC Voltammetry (SHACV)

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